

MINXRAY[®] INC.

DIGITAL IMAGING

HF120/60HPPWV POWER PLUS™

TECHNIQUE CHART for MinXray CR and DR IMAGING SYSTEMS*

VIEW	KvDC	mAs	Grid
Skull, AP	80	20	Y
Skull, Lat	80	12	Y
AP Chest	90	3	Y
Lat Chest	110	5	Y
Shoulder	70	4	N
Elbow, AP	64	2	N
Hand, AP	56	1	N
Hand, Lat	56	1	N
C-Spine	80	6	N
T-Spine, AP	90	12	Y
T-Spine, Lat	96	15	Y
L-Spine, AP	96	18	Y
L-Spine, Lat	96	22	Y
Abdomen, AP	80	12	Y
Pelvis, AP	96	18	Y
Hip, Lat	90	12	Y
Femur	80	10	N
Knee	70	3	N
Ankle	60	1.5	N
Foot, AP	60	1	N
Foot, Lat	64	1.35	N

SID = 40"

Grid is 6:1/178 lpi or equivalent

** While final results depend on many factors, if the SID is controlled per above and the grid used is as noted, the size of patients is the major reason for modifying a technique. The techniques noted are for an adult in the range of 160-185 lbs. If the patient is larger, increase the mAs; if smaller, decrease the mAs.*