

MinXray
HFFM (HFP40)
HFPS (HFP24)

Preliminary Podiatric TECHNIQUE CHART

Projection	CM	Distance	kV/mA	Time (Sec)
DP Foot 15° 0°	3 to 5	24"	60/10	0.04
	6 to 9	24"	60/10	0.05
	10 to 13	24"	60/10	0.06
Lat. Foot 90°	3 to 5	24"	60/10	0.06
	6 to 9	24"	60/10	0.07
	10 to 13	24"	60/10	0.08
Lat. Oblique 40°	3 to 5	24"	60/10	0.04
	6 to 9	24"	60/10	0.05
	10 to 13	24"	60/10	0.06
Med. Oblique 25°	3 to 5	24"	60/10	0.05
	6 to 9	24"	60/10	0.06
	10 to 13	24"	60/10	0.07
Lat. Hallux 90°	2	24"	60/10	0.04
	3 to 5	24"	60/10	0.05
	6 and above	24"	60/10	0.06
Axial Sesamoid 90°	Small	24"	60/10	0.08
	Medium	24"	60/10	0.09
	Large	24"	60/10	0.10
Axial Calcaneal 25°	4 to 6	24"	60/10	0.08
	7 to 10	24"	60/10	0.09
	11 to 14	24"	60/10	0.10
Harris & Beath 35° 45°	5 to 7	24"	60/10	0.10
	8 to 11	24"	60/10	0.11
	12 to 15	24"	60/10	0.12
AP/Mortice 90°	4 to 6	24"	60/10	0.07
	7 to 10	24"	60/10	0.08
	11 to 14	24"	60/10	0.09
Lateral 90°	3 to 5	24"	60/10	0.06
	6 to 9	24"	60/10	0.07
	10 to 13	24"	60/10	0.08
Lat. Oblique 90°	4 to 6	24"	60/10	0.07
	7 to 10	24"	60/10	0.08
	11 to 14	24"	60/10	0.09
Med. Oblique 90°	4 to 6	24"	60/10	0.07
	7 to 10	24"	60/10	0.08
	11 to 14	24"	60/10	0.09

Notes:

*Cassette/Film/Screen - Kodak Lanex Fine/100 R/S Speed Screens
 These are recommended starting techniques. Final results depend upon many factors. Therefore, if films are too dark (overexposed), decrease timer one stop; if films are too light (underexposed), increase timer one stop.*