

HF120/60HPPWV Power Plus™ /Orex ACL4 TECHNIQUE CHART

X-ray Unit: MinXray HF120/60HPPWV PowerPlus™/Orex ACL4

VIEW*	kVDC	mAs	400 SPEED SCREEN/FILM	GRID
			IMAGING SYSTEM TIME (sec.)	
AP Skull	76	5.92	0.18	Y
Lat Skull	70	3.15	0.07	Y
AP Chest	90	1.98	0.06	Y
Lat Chest	110	2.28	0.12	Y
Shoulder	70	2.25	0.05	N
Elbow	64	1.35	0.03	N
AP Hand	56	1.00	0.02	N
Lat Hand	56	1.00	0.02	N
C-Spine	80	3.04	0.08	Y
AP T-Spine	90	3.30	0.10	Y
Lat T-Spine	96	4.68	0.18	Y
AP L-Spine	96	4.47	0.17	Y
Lat L-Spine	96	9.30	0.40	Y
AP Abdomen	80	8.05	0.26	Y
AP Pelvis	96	4.47	0.17	Y
Lat Hip	90	3.84	0.14	Y
Femur	80	3.04	0.08	N
Knee	70	2.25	0.05	N
Ankle	60	1.00	0.02	N
AP Foot	60	1.00	0.02	N
Lat Foot	64	1.35	0.03	N

SID = 40"

Grid is 6:1/178 lpi

MINXRAY[®] INC.

3611 Commercial Avenue
Northbrook, Illinois 60062-1822, USA
1-847-564-0323
Toll Free 1-800-221-2245
FAX: 1-847-564-9040
E-mail: info@minxray.com
Web: www.minxray.com